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
Filter: Any Substrate ▼ Any Platform ▼

The table below shows data for the cell type and Nucleofector™ Platform selected. Those data are either based on Lonza Optimized Protocols  or on results shared from customers who performed an optimization based on our guidelines. In case no data are shown for the selected Nucleofector™ Platform, please take a look at our [optimization strategy](#) to get further guidance on how to easily determine optimal Nucleofection conditions yourself.

Protocol	Kit	Program	Cells	Efficiency	Viable Cells	Substrate	Format	Platform	
	T cell, human	T-023	5e6	36-58%		Plasmid (general)	1 µg	100 µl	I/II/2b
	P3	EO-115	1e6	63-77%	47-72%	Plasmid (general)	0.4 µg	20 µl	4D X-Unit
	P3	96-EO-115	1e6	63.3-77.3%	47.3-71.3%	Plasmid (general)	0.4 µg	20 µl	Shuttle

Lonza

T cell

T cell, human stim.

Peripheral blood lymphocytes; pre-stimulated.

Cell Type: T Cell
Tissue Origin: blood
Species: human

Cell Characteristics: Suspension

Transfection Information

 Lonza Optimized Protocol
 Optimization Guideline

Filter: Any Substrate ▼ Any Platform ▼

The table below shows data for the cell type and Nucleofector™ Platform selected. Those data are either based on Lonza Optimized

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- **Experiments** – display user-defined experiments
- **Custom programs** – create and manage custom Nucleofection™ Programs
- **Results** – reopen or load result files onto a USB stick plugged into the USB port in front of the Core Unit
- **Lonza programs** – version information about the cell type list
- **Display and audio** – adjust brightness of the touch screen display,time to activate standby and deactivate the touch tone
- **Date and time** – set and format date and time
- **Language and keyboard** – change language settings for display functions and keyboard
- ... more: To switch between pages of the settings list
- **Firmware update** – update the system software (see chapter 2.10 for details)
- **Device cleaning** – opens the tray of the functional modules to remove it for cleaning
- **System restart** – restart system
- **Transport lock** – activates the transportation lock of the system

-
- **Synchronize** – supports synchronization of the actual firmware with the PC Editor
 - **Version** – displays firmware version and serial number of all Modules

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